



Menus

Vigilens designer




This quick menu offers some common printing actions, as well as a Cursor icon to deselect the current tool.


Home



This menu offers :

- classic printing and editing tools
- default font setup
- report grid configuration






The grid is “magnetic”: when objects are moved with the mouse, they “stick” to the grid, whether it is displayed or not. If one wants to move an object finely without it being on the grid, it is necessary to use the options of the object (see below).


Body

RSORGN	RSORGN	SDDOCO	XHSHP
7778	7778	7778	7778

- Align
- Adjust
- Center within section
- Order
- Conditional format...
- Cut Ctrl+X
- Copy Ctrl+C
- Paste Ctrl+V
- Delete Del
- Properties...




Display



Display


You can have a cover page added to your report. These buttons toggle between editing the cover page and the data page.




the last button is used to display or not the grid.




Grid on:



Grid off:




Tool Bar




Hide or show [Sections panel](#) and [Database fields panel](#)

Zoom




Choice of zoom level, with the possibility of automatic adjustment to display either the entire page on the screen, or adapt the zoom level to have the full width of the report in the editing screen.

Full width zoom :



Full page zoom :




Format



Fields basic format parameters: font options, text alignment, conditional formatting.


To access all the parameters of a field, select it > right click > Properties or go to the Components menu > Properties

Component




[Dedicated page](#)

Layout



Align



These buttons can only be used when several items are selected simultaneously. To make a multiple selection, hold down the shift button and click successively on all the




objects of the selection. Only objects from one and the same section or sub-section can be selected simultaneously.



The **last** selected object is used as the reference.


Setup




Make width, height or both equals.

As for Align tool, this menu requires to have at least 2 objects selected, the last one being the reference.

Center within section



Order



From:
<https://vigilens.wiki/dokuwiki/> - Vigilens Reporting Knowledge Garden

Permanent link:
https://vigilens.wiki/dokuwiki/doku.php?id=en:v8_0_0:designer:menus:accueil:nouvellerequete:editeurimpression:menus:start

Last update: 2020/07/15 14:29

